

## Sentes

Ultra-efficient reclining bathing system, suitable for long-term care and acute care environments.



Designed to accommodate a diverse range of client needs, including clients with more challenging conditions, the Sentes bath offers care providers a versatile, ultra-efficient bathing solution.

The reclining bathing system and moulded seat help to improve postural control, enhancing bathing comfort and delivering the optimum blend of functionality and performance.

The vertical-raising, side-entry door also ensures easy-access for bathers, making transfers safer and simpler for care staff.

With its quick-fill technology and low water usage, it reduces time between refills, allowing carers to focus on the needs of bathers. Therefore, the Sentes is also extremely economical to run, particularly in high-traffic bathing areas.



### Key features:

- Accommodates a variety of different client needs with sculpted seat for improved postural control and reclining action
- Vertical-raising, side-entry door ensures easier access for bathers
- Ultra-efficient design with rapid-fill technology ensures low water usage to help reduce energy costs
- Bathers can transfer from standing or alternatively a ceiling track hoist or mobile hoist can be used, offering utmost flexibility in a multiuser application.
- Sleek, elegant bath design complements any bathroom
- Internal grab bar and head cushion ensure enhanced safety and comfort
- TMV3 controls maintain a safe, comfortable water temperature, reducing any risk of scalding
- Easy-to-use, one-touch tilt and recline operation
- Battery backup safety system in case of mains power failure
- Left handed bath configuration only
- Free standing design making it easy to install
- 205kg (32 stone) safe lifting weight accommodates the vast majority of bathers



# Sentes

## Key features and options



Easy to use recline operation



Lightweight, ergonomic hand operated shower unit option



Sculped seat design with hygiene aperture



Padded head cushion providing extra comfort



Easy to use controls



Chromotherapy lighting



### Key features and options:

- Hydrotherapy multi-speed air spa system
- KIWA compliant shower pole, (CAT 5) - Including shower head, diverter, spout and TMV3 valves
- Chromotherapy, multi-coloured LED lighting visually enhances the bathing experience
- Sound system – wirelessly controlled via Bluetooth from a smart device
- Autofill, thermostatically-controlled auto-fill system automatically fills the bath to a pre-set depth



# Sentes

## Technical specification

Length (mm)	1854
Width (mm)	865
Height (mm)	1295-1980
Footwell width (mm)	610
Bath height with shower pole (mm)	1536-1836
Maximum user weight	205kg
Water capacity	303 litres
Approximate bath weight (empty)	140kg
Approximate bath weight (full)	442kg
Typical mechanical operation times	38 secs

Class I Electrical Device	
Power supply	Mains 230V to be supplied via 30mA RCD
Operating voltage	24V
Battery backup	24V DC lead acid, for emergency use only
Actuator duty cycle	Max. 10% or 2min continuous use followed by 18min not in use.  The life cycle test has been performed with a stabilised power supply (10% duty cycle) on a 20mm stroke actuator at max. load for 10,000 cycles at ambient temperature
Water supply	Unblended Hot & Cold water at 1bar Min > 3.5bar Max pressure
IP Rating	IPX3
A' weighted sound emission	<60dB(A)
A' weighted sound emission with air spa system (option)	<75dB(A)
Ambient temperature range	5°C to 50°C
Relative humidity range	10-80%
Atmospheric pressure range	50kPa to 106kPa
Medical Device Category	Class I

### Antimicrobial protection

This bath includes BioCote® technology which has been proven to reduce bacteria and even some viruses, including influenza A H1N1, E.coli, MRSA and Pseudomonas.

### Intended use

This equipment is intended to facilitate the bathing of residents in care facilities or in a home care environment. The product is designed to help people who struggle with getting in and out of a bath, especially those with physical limitations. It's useful for both hygiene and relaxation purposes.

### Caution

To ensure best practice and safety the bath must be operated only by a trained professional with adequate understanding of the specific environment, its common practices and processes. It is vital that regular assessments of the bather are undertaken to ensure the bath remains suitable.

Please be aware of the sound emission and any optional accessories working attributes to ensure bather suitability before operation.

If you have any questions regarding this product or operation of the bath, feel free to contact your authorised representative or Gainsborough Healthcare Group Ltd. on 01527 400022.



**Gainsborough Healthcare Group**  
Life enhancing bathing.

01527 400022

info@gainsboroughbaths.com  
www.gainsboroughbaths.com

# Sentes

## Technical Drawings

### Sentes

Dimensions in millimetres (mm).

